

Technical Specifications

Dimensions and Weight

Description	Specification	
Injector (Pedestal)	Dimensions (WxDxH)	76.2 x 68.6 x 167.6 cm
	Weight (approximate)	76 kg
Control Room Unit	Dimensions (WxDxH)	29.2 x 20.2 x 25.9 cm
	Weight	7.2 kg

Injector Specifications

Description	Specification
Fluid Delivery Method	Piston-based
Fluid Loading	Contrast Media: 50 mL – 500 mL Saline: 50 mL – 1000 mL
Reservoirs (volume)	3 x 200 mL reservoirs (two same API Contrast Media and one saline)
Air Detection	3 inlet air sensors to detect when bottles and bags run empty 1 outlet air sensor to detect air before it reaches the Patient Line Reservoir air detection and automatic air removal
Battery Performance	Up to 16 hours or up to 60 injections
Displays	2 independent, Wi-Fi enabled, full color, touch screen units with modern graphical user interfaces located in the scan room and the control room
Barcode Reader	Integrated barcode reader on injector unit for contrast tracking and recording
Integrated Prime Container (volume)	680 ml

Protocols Specifications

Description	Specification
Volume Range	1 mL to 200 mL for contrast and saline phases; 1 mL to 400 mL for DualFlow phases (in 1 mL increments)
Flow Rate Range	0.1 to 10 mL/sec in 0.1 mL/sec increments
DualFlow	Simultaneous injection of contrast and saline in increments of 5%
Programmable Pressure Limit	From 50 psi to 300 psi in 1 psi increments (also displayable in kPa and kg/cm ²)
Programmable Delay	1 second to 900 seconds in 1 second increments
Pause	Maximum pause time of 20 minutes
Injection Capabilities	Maximum of 6 phases per injection Maximum of 10 injections per exam
Skip to Next Phase	Enables the user to stop the current phase of an injection and start the next phase
Patency Check / Test Inject	Real time flow rate adjustment and multi-phase test injections

Disposables

Description	Specification
Patient Line	Single use, 250 cm, snap-in single hand design, 2 check valves, auto-primers upon insertion, recessed connection to prevent accidental touch
Day Set	Pre-assembled, labeled for 24-hour use, designed for set-up in less than 2 minutes, recessed connection to prevent accidental touch
Replacement Spikes	Labeled for 24-hour use, can be changed in case of an accidental touch, without having to change the entire Day Set

Workflow Solutions

Description	Specification
Automated Documentation	For access to contrast and injection information in PACS: Integrated barcode reader, Modality Worklist, PACS Interface
Processing Productivity	RIS Interface, Speech Recognition Interface, Radimetrics Interface
Scanner Connectivity	ISI 700, ISI 900, Connect.CT
Smart Protocols	P3T algorithms calculate individualized injection protocols for cardiac CT, CT pulmonary angiography and CT imaging of abdominal organs Protocol Assistant provides basic protocol calculations for a range of procedure types
eGFR Calculator	Available on injector and Control Room Unit as part of the standard configuration

Bayer, the Bayer Cross, MEDRAD Centargo, MEDRAD, Centargo and P3T are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries. Other trademarks and company names mentioned herein are properties of their respective owners and are used herein solely for informational purposes. No relationship or endorsement should be inferred or implied.

© 2022 Bayer. This material may not be reproduced, displayed, modified or distributed without the express prior written consent of Bayer.



Bayer AB – Sverige
Berzelius väg 35
Box 606
SE-169 26 Solna
Tel: +46 (0) 31 748 28 80

Bayer A/S – Danmark
Arne Jacobsens Allé 13
DK-2300 København S
Tel: +45 38 16 16 16

Bayer AS – Norge
Drammensveien 288
Postboks 193
NO-1325 Lysaker
Tel: +47 22 06 57 10

More information on
radiology.bayer.se